

GC Derivatizing & Silylating Reagents

For Improved GC Separation

Over the years with the advancement in use of analytical instrument and data analysis platform like Gas Chromatography, GC-Derivatization has become an integral part of many Biotechnology, Diagnostic and Biopharmaceutical companies as a part of their sample preparation routine for GC analysis. SRL aid such industries by offering high quality GC-Derivatizing & Silylating Reagents, thereby supporting in better separation with increased resolution and response.

Prdt. Code	Product Name	CAS No.	Assay	Pkg. Unit
26477	N,O-Bis-(Trimethylsilyl) Acetamide (BSA)	10416-59-8	min. 95%	25ml, 100ml, 250ml
57677	N,O-Bis-(Trimethylsilyl) Trifluoro Acetamide (BSTFA)	25561-30-2	min. 99%	5ml, 25ml, 100ml
32431	tert-Butylchlorodiphenylsilane (TBDPSCI)	58479-61-1	min. 98%	10ml, 25ml, 100ml, 500ml
90611	N-tert-Butyldimethylsilyl Imidazole (TBDMSIM)	54925-64-3	min. 97%	1ml, 5ml
52573	tert-Butyldimethylsilyl Trifluoromethanesulfonate (TBDMS triflate)	69739-34-0	min. 98%	5ml, 25ml, 100ml
21800	Dimethyl Dichlorosilane (DMDCS)	75-78-5	min. 98%	500ml, 1000ml
28437	1,1,1,3,3,3-Hexamethyldisilazane (HMDS)	999-97-3	min. 98%	25ml, 100ml, 500ml, 2500ml
92009	N-Methyl-N-Trimethylsilyl Trifluoroacetamide (MSTFA)	24589-78-4	min. 98%	5ml, 25ml
13573	N-Methoxymethyl-N-Trimethylsilylmethyl-Benzylamine	93102-05-7	min. 98%	5gm, 25gm
27329	o-Phthalaldehyde (OPA)	643-79-8	min. 99%	5gm, 25gm, 100gm
49601	N-(3-Trimethoxysilylpropyl) Diethylenetriamine (TMSPD)	35141-30-1	min. 95%	100ml
20139	Trimethyl Chlorosilane (TMCS)	75-77-4	min. 98%	250ml, 1000ml
90265	Trimethyl Chlorosilane (TMCS)	75-77-4	min. 99%	100ml, 250ml, 1000ml
49901	N-(Trimethylsilyl)diethylamine (TMSDEA)	996-50-9	min. 97%	25ml, 100ml
45238	N-Trimethylsilylimidazole (TMSIM)	18156-74-6	min. 98%	5ml, 25ml, 100ml
97039	3-Trimethylsilylpropionic Acid (TMSPA)	5683-30-7	min. 95%	1ml, 5ml



Please contact:

SRL® SISCO
RESEARCH
LABORATORIES
PVT. LTD.

www.srlchem.com