



27839

## PreStained High Range 1 Tri-Coloured Protein Ladder: 10-180kDa

Part E

### Specifications

Appearance (Form)	Liquid
Appearance (Colour)	Dark Blue
Appearance (Clarity)	Clear
Electrophoresis test	Ladder should show 10 prestained protein bands during SDS PAG electrophoresis.

### Other Information

Description

BioLit Prestained High Range 1 Tri-Coloured Protein Ladder is a Ready-To-Use solution of carefully selected 10 stained proteins within a wide range of molecular mass (10-180 kDa). The ladder has 10 bands at 10,15,25,35,45,60,75,100,140 & 180 kDa covalently bonded with a chromophore. All the bands are bonded to a blue chromophore except for 75 kDa (red chromophore and 25 kDa (green chromophore) The three colour identification in ladder helps in easy orientation of the molecular sizes of the protein bands. The markers have been carefully prepared under sterile conditions and stored in appropriate additives to ensure long term stability. SDS PAGE analysis  
Western Blotting  
Allows monitoring protein migration during electrophoresis.

Application

### General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
0.25 ml	38229090 (GST 12%)
1 ml	38229090 (GST 12%)

## Available Packages

0.25 ml

1 ml

---

### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: [info@srchem.com](mailto:info@srchem.com), website [www.srchem.com](http://www.srchem.com)

---