

## 37351

# RNA Mini Spin Column Set (2ml SC w/ CT each)

Part L

			ica		
	$\mathbf{n}$	$\sim$ 1 $\pm$	-	+10	
_		$\sim$ 111			ш

Appearance (Form)

Absence of nucleases

Absorbance

Stability

Circular Columns and Collection tubes

Nucleic acids are incubated in the tube at 37°C for 1 hour

20-40µg

Stable if stored as recommended

#### **Other Information**

Description

A purification aide utilizing silica membrane in a column form is used for purifying/extracting nucleic acids of high quality. These columns are compatible with in-house buffers as well as store-brought ones.

The set is supplied with 2ml collection tubes.

### **General Information**

Storage	25 to 40°C (Room Temperature)
Shelf Life	36 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Nos	38229090 (GST 12%)
50 Nos	38229090 (GST 12%)

#### **Available Packages**

50 Nos

100 Nos