



50255

ProxiB 100bp DNA Ladder Plus

Part E

Specifications

Appearance (Form)	Liquid
Appearance (Colour)	Blue
A260/A280	1.7~19
Concentration	100-150ng/μl
Stability	No obvious degradation was seen after incubation at 37°C for 1 week. Well-defined bands are formed during agarose gel electrophoresis.

Other Information

Description

The ProxiB 100bp ladder plus is ideal for determining the size of double-stranded DNA from 100 to 3000 base pairs. The ladder consists of 14 linear double-stranded fragments. The 500bp and 1200bp fragment are present at increased intensity to allow easy identification. All fragments are precisely quantified and mixed during the production. For 5 mcl loading, all fragments except 500bp and 1200bp are 40ng. The 500bp and 1200bp fragment are 100ng. This ladder is pre-mixed with loading dye and is ready to use.

Includes

50 mcg
• ProxiB 100bp Ladder Plus - 50 mcg

5 x50 mcg
• ProxiB 100bp Ladder Plus - 5 x 50 mcg

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 x50 mcg	38229090 (GST 12%)
50 mcg	38229090 (GST 12%)

Available Packages

50 mcg

5 x50 mcg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
