



54896

## FluroBronze 500bp Ladder

Part E

### Specifications

Appearance (Form)	Liquid
Appearance (Colour)	Blue
A260/A280	1.7~1.9
Concentration	~90ng/μl
Stability	No obvious degradation was seen after incubation at 37°C for 1 week. Well-defined bands are formed during agarose gel electrophoresis.

### Other Information

Description

The FluroBronze 500bp ladder is specially prepared as ready-to-load fluorescent ladders ideal for determining the size of double-stranded DNA from 500 to 5,000 base pairs. The ladder consists of 10 linear double-stranded fragments. The 2,500bp fragment is present at increased intensity to allow easy identification.. The flurobronze ladder contains in addition to the tracking dye, a non-carcinogenic nucleic acid stain. The ladder should be loaded onto an agarose gel (of appropriate concentration) containing no fluorescent dye. The stain contained in the ladder replaces ethidium bromide, thus eliminating its addition to the gel.

- FluroBronze 500bp Ladder - 50 mcg
- FluroBronze Stain - 250 mcl

Includes

### General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
50 mcg	38229090 (GST 12%)

### Available Packages

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com