



68738

3,3,5,5-Tetramethylbenzidine (TMB) extrapure
AR, ExiPlus, Multi-Compendial, 98%

Part B

CAS : 54827-17-7

Molecular Formula : C₁₆H₂₀N₂

Molecular Weight : 240.35

Specifications

ExiPlus	meets compendial specs of BP, Ph.Eur
Appearance (Colour)	Off - white to brown
Appearance (Form)	Crystalline powder
Assay (NT)	min. 98.0%
Melting Point	168 - 170°C

Other Information

3,3,5,5-Tetramethylbenzidine (TMB) is a colorimetric substrate used with peroxidase and peroxidase-coupled systems. It produces a soluble end product, which is blue in color and can be read spectrophotometrically at 370 nm or 655 nm(photometric mode). The reaction can be stopped with acid (Trichloroacetic acid or Sulphuric Acid or as per application). This product is also been used for hematology and immunochemical tests.

Preparation Instructions:

A TMB stock solution can be prepared by dissolving 10 mg/ml in DMSO(41.6 mM sol). The working solution is made by diluting 0.1 ml of the DMSO stock with 9.9 ml Sodium Citrate buffer, pH 6.0,(sodium Acetate (0.7M) and adjust ph by using citric Acid 6.0) and then adding 2ul of fresh 30% hydrogen peroxide per 10 ml of substrate/buffer solution immediately prior to use. The hydrochloride salt can be directly dissolved in aqueous media. For this reason, the hydrochloride salt is the preferred form of TMB; however, both forms will give identical results.

Stock solutions of TMB in DMSO can be frozen in 1–2 ml aliquots at –20deg C. Solutions may appear to remain active for at least one years when stored in this manner. Any oxidation would be revealed as a color change in the working solution.

Application

General Information

Storage2 to 8°C (Refrigerate)

Shelf Life60 Months

IMDG IdentificationNot Regulated for Transport (Non-Haz)

HSN Code

1 Gms29215990 (GST 18%)

Available Packages

1 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com