

83141

Bacteroides Selective Supplement (50mg/vl)

Part D

Specifications
Appearance (Colour)
Solubility
Cultural Response
Organism
Inoculum (cfu) 10-100
Organism
Inoculum (cfu) 10-100
Organism
Inoculum
Organism
Inoculum

Other Information

Applications

Used for the selective isolation of Bacteroides species.

Composition

Gentamicin 50.00mg/vial

Directions

- 1. Aseptically add 5ml of sterile distilled water to 1 vial.
- 2. Mix well and aseptically add to sterile cooled media at 45-50°C and pour in sterile petriplates.
- 3. Add 5ml (1 vial) of supplement in 500 ml of Bacteroides Bile Esculin Agar Base.

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
5 vl	38229090 (GST 12%)

Available Packages 5 vl

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.