



52075

ProxiO 100bp DNA Ladder Plus

Part E

Specifications

| | |
|---------------------|--|
| Appearance (Form) | Liquid |
| Appearance (Colour) | Blue |
| A260/A280 | 1.7~1.9 |
| Concentration | 120 - 150ng/ μ l |
| Stability | No obvious degradation was seen after incubation at 37°C for 1 week. Well-defined bands are formed during agarose gel electrophoresis. |

Other Information

Description

The ProxiO 100bp ladder plus is ideal for determining the size of double-stranded DNA from 100 to 3000 base pairs. The ladder consists of 14 linear double-stranded fragments. The 500bp and 1200bp fragment are present at increased intensity to allow easy identification. All fragments are precisely quantified and mixed during the production. For 5 mcl loading, all fragments except 500bp and 1200bp are 40ng. The 500bp and 1200bp fragment are 100ng

Includes

50 mcg
• ProxiO 100bp Ladder Plus - 50 mcg
• 6X Orange Loading Dye - 1.0 ml

5 x50 mcg
• ProxiO 100bp Ladder Plus - 5 x 50 mcg
• 6X Orange Loading Dye - 5 x 1.0 ml

General Information

| | |
|---------------------|---------------------------------------|
| Storage | -20 °C (Blue/Dry Ice) |
| Shelf Life | 24 Months |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code | |
| 5 x50 mcg | 38229090 (GST 12%) |
| 50 mcg | 38229090 (GST 12%) |

Available Packages

50 mcg

5 x50 mcg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
