



10253

3,5-Dinitrosalicylic Acid pure, 98%

Part B

CAS : 609-99-4

Molecular Formula : C7H4N2O7

Molecular Weight : 228.19

Specifications

Appearance (Colour)	Pale yellow
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% solution in methanol	Clear
Solubility (Colour) 10% solution in methanol	Yellow
Assay (NaOH titration)	min. 98%
Melting Point	169 - 175°C

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
HSN Code	
25 Gms	29182190 (GST 18%)
100 Gms	29182190 (GST 18%)
500 Gms	29182190 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle

Available Packages

25 Gms
100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com