



10344

## Low-Range DNA Marker: 0.1 – 1 kb

Part E

### Specifications

Appearance (Form)	Liquid
Appearance (Colour)	Blue
Appearance (Clarity)	Clear
Electrophoresis test	Well-defined bands are formed during agarose gel electrophoresis
Stability	No obvious degradation was seen after incubation at 37°C for 1 week

### Other Information

Description

This marker range has been designed to assist in determining the size of PCR product, insert and restriction analysis. This marker contains Bromophenol Blue as the tracking dye.

### General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 mcg	38229090 (GST 12%)
15 mcg	38229090 (GST 12%)
50 mcg	38229090 (GST 12%)

### Available Packages

15 mcg  
50 mcg  
100 mcg

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com