



11249

L-Arginine Methyl Ester Dihydrochloride extrapure, 98%

Part B

CAS : 26340-89-6

Molecular Formula : C₇H₁₆N₄O₂.2HCl

Molecular Weight : 261.15

Specifications

Appearance (Colour)	White to pale yellow
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Assay	min. 98%
Melting Point (Decomposition)	~190°C

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29379020 (GST 18%)
25 Gms	29379020 (GST 18%)
5 Gms	29379020 (GST 18%)
Type of Packing	
25 Gms	Glass Bottle

Available Packages

5 Gms
25 Gms
100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com