



11287

Morpholine extrapure AR, ACS, ExiPlus, Multi-Compendial, 99.5%

Part A

CAS : 110-91-8

Molecular Formula : C₄H₉NO

Molecular Weight : 87.12

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Colour (APHA)	max. 15
Assay	min. 99.5%
Density (g/ml) @ 20°C	0.999-1.000
Refractive Index (20°C)	1.454-1.455
Boiling Range	127-129°C
Non Volatile Matter	max. 0.01%
Iron (Fe)	max. 0.00005%
Heavy Metals (Pb)	max. 0.0001%
Copper (Cu)	max. 0.00001%
Water (KF)	max. 0.1%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2054 , IMCO Class No.:8(3) , Packing Group:I
HSN Code	
500 ml	29349920 (GST 18%)
Type of Packing	
500 ml	Glass Bottle

Available Packages

500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com