



11443

Dish Evaporating Porcelain, 20ml, OD 50mm, RB, Deep with Spout

Part L

Specifications

Appearance (Colour)	White to Off white
Appearance (Form)	Evaporating dish
Material of Construction (MOC)	Porcelain
Evaporating Dish Design Type	Deep with Spout, Round bottom
Evaporating Dish Top Outer Dia	50 mm \pm 2.5mm
Evaporating Dish Height	20.5mm \pm 5mm
Evaporating Dish Weight	21.5gm \pm 5gm
Max Usable Temperature	1150°C \pm 50°C
Flextural Strength @20°C	Min. 80 Mpa (815 Kg/cm ²)

Other Information

Application	A high quality evaporating basin for critical laboratory use.
-------------	---

General Information

Storage	25 to 40°C (Room Temperature)
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
6 Nos	69091100 (GST 18%)

Available Packages

6 Nos

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com