



11814

Elastase (ELA, PPE) ex. Porcine Pancreas,  
10U/mg Protein

Part B

CAS : 39445-21-1

Specifications

Appearance (Colour)	White to off white
Appearance (Form)	Powder
Activity	min 10 units/mg protein
Protein content	55 to 85%
Unit Definition	One unit of Elastase will hydrolyze 1.0 µmole of N-succinyl-L-Ala-Ala-Ala-p-nitroanilide per minute at pH 8.0 at 25 °C.
Enzyme solution for activity determination	Prepared in 100 mM Tris HCl, pH 8.0.

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Mg	35079099 (GST 18%)

Available Packages

10 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com