



11996

Sodium Dichromate Dihydrate extrapure AR, ACS, ExiPlus, Multi-Compendial, 99.5%

Part A

CAS : 7789-12-0

Molecular Formula : Na₂Cr₂O₇.2H₂O

Molecular Weight : 298.00

Specifications

ExiPlus	meets compendial specs of USP
Appearance (Colour)	Orange
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Orange
Assay (Iodometric)	min. 99.5%
Insoluble matter	max. 0.005%
Chloride (Cl)	max. 0.001%
Sulphate (SO ₄)	max. 0.01%
Iron (Fe)	max. 0.002%
Heavy Metals (Pb)	max. 0.002%
Copper (Cu)	max. 0.001%
Potassium (K)	max. 0.01%
Calcium (Ca)	max. 0.003%
Aluminium (Al)	max. 0.002%
Magnesium (Mg)	max.0.005%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3086 , IMCO Class No.:6.1(5.1) , Packing Group:I
HSN Code	
5 Kg	28413000 (GST 18%)
500 Gms	28413000 (GST 18%)
Type of Packing	
500 Gms	Plastic Bottle

Available Packages

500 Gms

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
