



12358

## p-Toluidine pure, 98%

Part B

CAS : 106-49-0

Molecular Formula : C<sub>7</sub>H<sub>9</sub>N

Molecular Weight : 107.16

### Specifications

Appearance (Colour)	Slight brownish
Appearance (Form)	Crystalline compound
Solubility (Turbidity) 10% solution in methanol	Clear
Assay (NT)	min 98%
Melting Point	42 - 44°C

### General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	UN No.:3451 , IMCO Class No.:6.1 , Packing Group:II
HSN Code	
500 Gms	29215990 (GST 18%)
Type of Packing	
500 Gms	Glass Bottle

### Available Packages

500 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com