



12375

## Agarose TBE mix 2.0%

Part B

### Specifications

Appearance (Colour)	White to off white
Appearance (Form)	Powder
Function Test?	Passes

### Other Information

Notes	Pre-mixed blend of Agarose and TBE at 1x concentration.
-------	---

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	38229090 (GST 12%)
Type of Packing	
100 Gms	Plastic Bottle

### Available Packages

100 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com