



12447

## Mid-Range DNA Marker: 0.5 – 5 kb

Part E

### Specifications

Appearance (Form)	Liquid
Appearance (Colour)	Blue
Appearance (Clarity)	Clear
Electrophoresis test	Well-defined bands are formed during agarose gel electrophoresis.
Stability	No obvious degradation was seen after incubation at 37°C for 1 week

### Other Information

Description

This marker range should facilitate the analysis of inserts cloned into commonly used plasmids and restriction digests. This marker contains Bromophenol Blue as the tracking dye.

### General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	24 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 mcg	38229090 (GST 12%)
15 mcg	38229090 (GST 12%)
50 mcg	38229090 (GST 12%)

### Available Packages

15 mcg

---

50 mcg

---

100 mcg

---

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com