



12475

Triethyl Citrate extrapure, 99%

Part A

CAS : 77-93-0

Molecular Formula : C₁₂H₂₀O₇

Molecular Weight : 276.29

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to slight yellow
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	1.136 - 1.138
Refractive Index (20°C)	1.442 - 1.444
Acidity (C ₆ H ₈ O ₇)	max. 0.01%
Water (KF)	max. 0.2%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	36 Months
HSN Code	
500 ml	29181590 (GST 18%)
1000 ml	29181590 (GST 18%)
2500 ml	29181590 (GST 18%)

Available Packages

250 ml

1000 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com