



12502

Polysorbate 60 (Tween 60) extrapure

Part A

CAS : 9005-67-8

Specifications

Appearance (Colour)	Pale yellow
Appearance (Clarity)	Slightly hazy
Appearance (Form)	Viscous liquid
Density (g/ml) @ 20°C	1.079 - 1.081g/ml
Acid value	max. 2.0
Hydroxyl value	81 - 96
Saponification value	45 - 55
Water (KF)	max. 2%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	34024900 (GST 18%)
5 Kg	34024900 (GST 18%)
500 Gms	34024900 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
500 Gms	Glass Bottle

Available Packages

100 Gms

500 Gms

5 Kg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com