



12590

20X SSC Buffer pH 6.9-7.1 for molecular biology

Part B

Specifications

| | |
|------------------------|--------------|
| Appearance (Colour) | Clear |
| Appearance (Clarity) | Colourless |
| Appearance (Form) | Solution |
| pH | 6.9 - 7.1 |
| DNase, RNase, Protease | Not detected |

General Information

| | |
|------------|------------------------|
| Storage | 2 to 8°C (Refrigerate) |
| Shelf Life | 36 Months |
| HSN Code | |
| 200 ml | 38229090 (GST 12%) |
| 500 ml | 38229090 (GST 12%) |

Available Packages

200 ml
500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com