



12644

## Pyridine-3-Sulphonic Acid (P3SA) extrapure, 98%

Part B

CAS : 636-73-7

Molecular Formula : C<sub>5</sub>H<sub>5</sub>NO<sub>3</sub>S

Molecular Weight : 159.16

### Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 1% aq. solution	Clear
Solubility (Colour) 1% aq. solution	Colourless
Assay (T)	min. 98%
Water (KF)	max. 2%

### General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
HSN Code	
25 Gms	29333990 (GST 18%)

### Available Packages

25 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com