



12680

Sodium Cyanoborohydride pure, 95%

Part A

CAS : 25895-60-7

Molecular Formula : NaBH₃CN

Molecular Weight : 62.84

Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Powder
Assay (NT) (On dry basis)	min. 95%
Loss on drying	max. 6%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:3134 , IMCO Class No.:4.3(6.1) , Packing Group:II
HSN Code	
100 Gms	29420090 (GST 18%)
25 Gms	29420090 (GST 18%)
5 Gms	29420090 (GST 18%)
Type of Packing	
100 Gms	Glass Bottle
25 Gms	Glass Bottle

Available Packages

5 Gms

25 Gms

100 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com