



12692

Phenol Tris Equilibrated for molecular biology w/ Stabilizer

Part A

CAS : 108-95-2

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Yellow
Appearance (Form)	Liquid
pH	7.8 - 8.2
DNase, RNase, Protease	Not detected
Stabilizer (8-Hydroxyquinoline)	Stabilized with 0.05% 8-Hydroxyquinoline

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	UN No.:2821 , IMCO Class No.:6.1 , Packing Group:II
HSN Code	
100 ml	29071190 (GST 18%)
500 ml	29071190 (GST 18%)
60 ml	29071190 (GST 18%)
Type of Packing	
60 ml	Glass Bottle

Available Packages

60 ml
100 ml
500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com