



13057

Crucible Porcelain, 50ml, Wide with Lid

Part L

Specifications	
Appearance (Colour)	White to off White
Appearance (Form)	Crucible
Material of Construction (MOC)	Porcelain
Crucible Dia	Top OD 61mm ± 3
Crucible Outer Height	37mm ± 5
Crucible Wall Thickness	1.8mm ± 0.3
Crucible Weight	35gm ± 5
Max Usable Temperature	1150°C ± 50
Flextural Strength @20°C	Min. 80 Mpa (815 Kg/cm2)
Lid Dia	ID 69mm ± 3

General Information	
Storage	25 to 40°C (Room Temperature)
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Nos	69091100 (GST 18%)
4 Nos	69091100 (GST 18%)

Available Packages
1 Nos
4 Nos

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com