



13492

p-Arbutin extrapure, 98%

Part B

CAS : 497-76-7

Molecular Formula : C₁₂H₁₆O₇

Molecular Weight : 272.26

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Colour) 10% aq. solution	Clear
Solubility (Turbidity) 10% aq. solution	Colourless
pH (1% aq. solution)	5.0 - 7.0
Assay	min. 98%
Melting Point	195 - 198°C
Specific rotation [?]	-63° to -66° (C=2, Water)

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Gms	35040099 (GST 18%)
25 Gms	35040099 (GST 18%)
50 Gms	35040099 (GST 18%)
Type of Packing	
25 Gms	Plastic Bottle

Available Packages

10 Gms

25 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com