



13627

Isopropenyl Acetate extrapure, 99%

Part A

CAS : 108-22-5

Molecular Formula : C₅H₈O₂

Molecular Weight : 100.12

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	0.919-0.920
Refractive Index (20°C)	1.400-1.401
Boiling Range	95-97°C
Acidity (CH ₃ COOH)	max. 0.1%
Water (KF)	max. 0.1%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2403 , IMCO Class No.:3 , Packing Group:II
HSN Code	
100 ml	29153999 (GST 18%)
500 ml	29153999 (GST 18%)

Available Packages

100 ml

500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com