



13674

Sulpholane Anhydrous extrapure AR, 99%

Part B

CAS : 126-33-0

Molecular Formula : C₄H₈O₂S

Molecular Weight : 120.17

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to almost colourless
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 25°C	1.261-1.270
Refractive index (25°C)	1.484-1.485
Melting Point	27-28°C
Water (KF)	max. 0.2%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
2500 ml	29309099 (GST 18%)
500 ml	29309099 (GST 18%)

Available Packages

500 ml

2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com