



13814

Xylanase ex. Aspergillus Niger, 60MU/g

Part B

CAS : 37278-89-0

Specifications

|  |   |
|--|---|
| Appearance (Colour)                        | Off-white to pale brown   |
| Appearance (Form)                          | Powder  |
| Activity                                   | min. 60 mu/g powder   |
| Unit Definition                            | One unit will liberate 1.0 micromole of reducing sugar measured as xylose equivalent from xylan per minute at pH 4.5 at 30°C. |
| Enzyme solution for activity determination | Prepared in citrate buffer, pH 4.5.   |

General Information

|                     |                                       |
|---------------------|---------------------------------------|
| Storage             | 0 to 4°C (Freeze)                     |
| Shelf Life          | 60 Months                             |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code            |                                       |
| 100 Gms             | 35079099 (GST 18%)                    |
| 25 Gms              | 35079099 (GST 18%)                    |
| 500 Gms             | 35079099 (GST 18%)                    |
| Type of Packing     |                                       |
| 100 Gms             | Plastic Bottle                        |

Available Packages

|         |
|---------|
| 25 Gms  |
| 100 Gms |
| 500 Gms |

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com