



14420

Dextran Blue 500 ex. Leuconostoc Sp.

Part B

CAS : 87915-38-6

Specifications

Appearance (Colour)	Blue
Appearance (Form)	Compound
Average Molecular weight	made from dextran 500 KDa
Loss on drying	max. 10%
Degree of substitution (mmol/g)	0.01 - 0.08

General Information

Storage	-20 °C (Blue/Dry Ice)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Gms	35051090 (GST 18%)

Available Packages

1 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com