



14443

Multipurpose Test Kit (Pre-Weighed Media)

Part D

Other Information

Applications

Used for isolation, identification and cultivation of a wide variety of microorganism.

Composition

(89912)Violet Red Bile Glucose Agar, Harmonized	41.53 gm/lt.
(50222)Baird Parker Agar(U/P),(I/P),(I)	63.00 gm/lt.
(65774)Brilliant Green Agar	55.09 gm/lt.
(21493)Columbia Agar Base, Harmonized	44.00 gm/lt.
90465 Cetrimide Agar Base, Harmonized	45.30 gm/lt.
(46675) EMB Agar Levine	37.46 gm/lt.
(37775)MacConkey Agar, Harmonized	50.03 gm/lt.

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	36 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
1 Kit	38229090 (GST 12%)

Available Packages

1 Kit

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com