



14749

Polyuridylic Acid Potassium Salt (POLY U Potassium) extrapure

Part B

CAS : 27416-86-0

Specifications

Appearance (Colour)	White
Appearance (Form)	Lyophilized powder
Solubility (Turbidity) 1% aq. solution	Clear
Solubility (Colour) 1% aq. solution	Colourless
Absorption maxima	258 - 262 nm
Spectral Ratio (A250nm/A260nm)	0.72 - 0.82
Spectral Ratio (A280nm/A260nm)	0.29 - 0.39
Spectral Ratio (A290nm/A260nm)	0.02 - 0.07
Water (KF)	Max. 10%

General Information

Storage	-20 °C (Blue/Dry Ice)
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
25 Mg	29349990 (GST 18%)

Available Packages

25 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com