



14887

Propidium Iodide Solution (1mg/ml in H₂O)

Part B

CAS : 25535-16-4

Specifications

| | |
|----------------------------|--------------------------|
| Appearance (Clarity) | Clear |
| Appearance (Colour) | Orange to reddish orange |
| Appearance (Form) | Liquid |
| Suitability for Microscopy | Passes test |

Other Information

Application

Suitable for use in fluorescence microscopy, laser scanning microscopy, flow cytometry, and fluorometry.

General Information

| | |
|---------------------|---------------------------------------|
| Storage | 2 to 8°C (Refrigerate) |
| Shelf Life | 24 Months |
| IMDG Identification | Not Regulated for Transport (Non-Haz) |
| HSN Code | |
| 10 ml | 38229090 (GST 12%) |
| Type of Packing | |
| 10 ml | Plastic Bottle |

Available Packages

10 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com