



15127

Phenyl Isothiocyanate (PITC) for HPLC & Sequencing, 99%

Part B

CAS : 103-72-0

Molecular Formula : C₇H₅NS

Molecular Weight : 135.19

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Almost colourless to yellow
Appearance (Form)	Liquid
Colour (APHA)	max. 50
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	1.130 - 1.132
Refractive Index (20°C)	1.648 - 1.650

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
HSN Code	
5 x1 ml	29309099 (GST 18%)
5 ml	29309099 (GST 18%)
5 x5 ml	29309099 (GST 18%)

Available Packages

5 ml
5 x1 ml
5 x5 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com