



15354

Sulphanilic Acid pure, 99%

Part A

CAS : 121-57-3

Molecular Formula : C₆H₇SO₃N

Molecular Weight : 173.19

Specifications

Appearance (Colour)	White to off - white to beige
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% solution in 5% Na ₂ CO ₃	Clear to opalescent
Assay (NaOH titration)	min. 99%
Loss on drying at 105°C	max. 0.5%
Iron (Fe)	max. 0.05%
Chloride (Cl)	max. 0.01%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29214234 (GST 18%)
500 Gms	29214234 (GST 18%)

Available Packages

100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com