



15360

Perovskite Quantum Dots (pQD 480),
Fluorescence 480nm (Hydrophobic)

Part C

Specifications

Appearance (Colour)	Green
Appearance (Form)	Crystals
Emission wavelength	480nm

Other Information

Notes	Readily soluble in toluene, chloroform and similar solvents. Not soluble in alcohols, ethers and polar solvents.
-------	--

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	36 Months
HSN Code	
50 Mg	28419000 (GST 18%)

Available Packages

50 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com