



15377

CdSe/ZnS Core Shell Quantum Dots, Fluorescence 630nm, (Hydrophobic)

Part C

Specifications

Appearance (Colour)	Brown
Appearance (Form)	Compound
Emission wavelength	630nm

Other Information

Notes	Hydrophobic, Highly luminescent semiconductor.
-------	--

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	36 Months
IMDG Identification	UN No.:3077 , IMCO Class No.:9 , Packing Group:III
HSN Code	
10 Mg	32065000 (GST 18%)
Type of Packing	
10 Mg	Plastic Vial

Available Packages

10 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com