



16024

Crucible Porcelain, 2ml, Wide with Lid

Part L

Specifications

Appearance (Colour)	White to off white
Appearance (Form)	Crucible
Material of Construction (MOC)	Porcelain
Crucible Dia	Top OD 18mm \pm 0.5
Crucible Outer Height	12mm \pm 0.5
Crucible Wall Thickness	1.3mm \pm 0.1
Crucible Weight	3gm \pm 0.1
Max Usable Temperature	1150°C \pm 50
Flextural Strength @20°C	Min. 80 Mpa (815 Kg/cm ²)
Lid Dia	ID 26mm \pm 2.6

General Information

Storage	25 to 40°C (Room Temperature)
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
6 Nos	69091100 (GST 18%)

Available Packages

6 Nos

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com