



16261

Disodium Tartrate Dihydrate extrapure AR, ACS, ExiPlus, Multi-Compendial, 99%

Part A

CAS : 6106-24-7

Molecular Formula : C₄H₄Na₂O₆.2H₂O

Molecular Weight : 230.08

Specifications

ExiPlus	meets compendial specs of USP, BP, Ph.Eur
Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 10% aq. solution	Clear
Solubility (Colour) 10% aq. solution	Colourless
Assay (NT)	min. 99%
pH (5% aq. solution)	7.0 - 9.0
Insoluble matter	max. 0.005%
Water (KF)	15.61 - 15.71%
Chloride (Cl)	max. 0.0005%
Sulphate (SO ₄)	max. 0.005%
Iron (Fe)	max. 0.0005%
Heavy Metals (Pb)	max. 0.0005%
Phosphate (PO ₄)	max. 0.0005%
Copper (Cu)	max. 0.0005%
Nitrogen Compounds (N)	max. 0.002%
Potassium (K)	max. 0.002%
Ammonium (NH ₄)	max. 0.003%
Calcium (Ca)	max. 0.01%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
250 Gms	29181390 (GST 18%)
500 Gms	29181390 (GST 18%)
Type of Packing	
250 Gms	Plastic Bottle

500 Gms

Plastic Bottle

Available Packages

250 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
