



16666

Isopropanol (IPA) Dried, 99.5%, water 0.005%

Part A

CAS : 67-63-0

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC)	min. 99.5%
Density (g/ml) @ 20°C	0.784 - 0.786
Refractive Index (20°C)	1.377 - 1.378
Boiling Range	81 - 83°C
Non Volatile Matter	max. 0.005%
Acidity (C3H6O2)	max. 0.001%
Water (KF)	max. 0.005%

General Information

Storage	25 to 40°C (Room Temperature)
IMDG Identification	UN No.:1219 , IMCO Class No.:3 , Packing Group:II
HSN Code	
500 ml	29051220 (GST 18%)
1000 ml	29051220 (GST 18%)

Available Packages

500 ml

1000 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com