



16881

## Z-L-Phenylalanine-4-Nitrophenyl Ester extrapure, 98%

Part B

CAS : 2578-84-9

Molecular Formula : C<sub>23</sub>H<sub>20</sub>N<sub>2</sub>O<sub>6</sub>

Molecular Weight : 420.41

### Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Amorphous powder
Solubility (Turbidity) 10% solution in DMF	Clear
Solubility (Colour) 10% solution in DMF	Colourless to pale yellow
Assay (HPLC) (on anhydrous basis)	min. 98 %
Specific rotation [?] @25°C	-24° to -26° (C=1,DMF)
Water (KF)	max.2%

### General Information

Storage	0 to 4°C (Freeze)
Shelf Life	60 Months
HSN Code	
10 Gms	29224990 (GST 18%)
25 Gms	29224990 (GST 18%)

### Available Packages

1 Gms  
10 Gms

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com