



17017

L-Serine Methyl Ester Hydrochloride extrapure, 99%

Part B

CAS : 5680-80-8

Molecular Formula : C₄H₉NO₃.HCl

Molecular Weight : 155.50

Specifications

Appearance (Colour)	White to off-white
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 5% aq. solution	Clear
Assay	min. 99%
Melting Point	160 - 165°C
Specific rotation [?]	+5.0 to +5.5° (C=2 , Methanol)

General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29379020 (GST 18%)
25 Gms	29379020 (GST 18%)
5 Gms	29379020 (GST 18%)
Type of Packing	
25 Gms	Glass Bottle
5 Gms	Glass Bottle

Available Packages

5 Gms
25 Gms
100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com