



17083

Phenyl Isothiocyanate (PITC) pure, 98%

Part B

CAS : 103-72-0

Molecular Formula : C₇H₅NS

Molecular Weight : 135.19

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Brownish yellow
Appearance (Form)	Liquid
Assay (GC)	min. 98%
Density (g/ml) @ 20°C	1.130 - 1.132
Refractive Index (20°C)	1.649 - 1.650

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2927 , IMCO Class No.:6.1(8) , Packing Group:II
HSN Code	
100 ml	29309099 (GST 18%)
500 ml	29309099 (GST 18%)
Type of Packing	
500 ml	Glass Bottle

Available Packages

25 ml

100 ml

500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com