



17798

Biuret extrapure AR, 97%

Part A

CAS : 108-19-0

Molecular Formula : CH₅O₂N₃

Molecular Weight : 103.08

Specifications

Appearance (Colour)	White to off white
Appearance (Form)	Crystalline powder
Assay (ex N)	min. 97% (on dry basis)
Water (KF)	max. 2 %
Melting Point (Decomposition)	184-190°C

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	29241900 (GST 18%)
25 Gms	29241900 (GST 18%)

Available Packages

25 Gms
100 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com