



17996

Aliquat 336, 88%

Part A

CAS : 63393-96-4

Molecular Formula : C25H54ClN

Molecular Weight : 404.16

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless to pale yellow
Appearance (Form)	Viscous liquid
Assay	min. 88%
Density (g/ml) @ 20°C	0.887 - 0.890 g/ml
Refractive Index (20°C)	1.466 - 1.467
Water (KF)	max. 5.0%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:2810 , IMCO Class No.:6.1 , Packing Group:III
HSN Code	
1000 ml	29239000 (GST 18%)
250 ml	29239000 (GST 18%)

Available Packages

250 ml
1000 ml
2500 ml
25 Ltr

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com