



18473

Bismuth (III) Sulphate pure, 90%

Part A

CAS : 7787-68-0

Molecular Formula : $\text{Bi}_2(\text{SO}_4)_3$

Molecular Weight : 706.15

Specifications

Appearance (Colour)	White
Appearance (Form)	Crystalline powder
Solubility (Turbidity) 5% solution in 5N HNO_3	Clear
Solubility (Colour) 5% solution in 5N HNO_3	Colourless
Assay (Complexometric)	min. 90%
Chloride (Cl)	max. 0.01%
Sodium (Na)	max. 0.01%
Potassium (K)	max. 0.005%

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
100 Gms	28332990 (GST 18%)
500 Gms	28332990 (GST 18%)
Type of Packing	
100 Gms	Plastic Bottle

Available Packages

100 Gms
500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srlchem.com, website www.srlchem.com