



18824

Chitosan (Medium MW) extrapure, 150-500m.Pas, 90% DA

Part B

CAS : 9012-76-4

Specifications

Appearance (Colour)	White to pale yellow
Appearance (Form)	Powder
Solubility (Colour) 1% solution in 1% Acetic acid	Colourless to pale yellow
Loss on drying	max. 5%
Viscosity (20°C)	150 - 500 m.Pas
Deacetylated Degree	min. 90%

Other Information

Application

Biocompatible, antibacterial and environmentally friendly polyelectrolyte with a variety of applications including water treatment, chromatography, additives for cosmetics, textile treatment for antimicrobial activity, novel fibers for textiles, photographic papers, biodegradable films, biomedical devices, and micro-capsule implants for controlled release in drug delivery.

General Information

Storage	8 to 25°C (Cool & Dry Area)
Shelf Life	60 Months
IMDG Identification	Not Regulated for Transport (Non-Haz)
HSN Code	
10 Gms	29211990 (GST 18%)
100 Gms	29211990 (GST 18%)
25 Gms	29211990 (GST 18%)
500 Gms	29211990 (GST 18%)
Type of Packing	
10 Gms	Glass Bottle
100 Gms	Glass Bottle
25 Gms	Glass Bottle
500 Gms	Plastic Bottle

Available Packages

10 Gms

25 Gms

100 Gms

500 Gms

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com
