



18989

SWCNT Type 4 Carbon Nanotubes Single-Walled, Length 5-30?m, OD 1-2nm, 90%

Part C

CAS : 308068-56-6

Specifications

| | |
|---------------------|----------|
| Appearance (Colour) | Black |
| Appearance (Form) | Powder |
| Assay | min. 90% |
| OD | 1 - 2nm |
| Length | 5 - 30?m |
| -COOH content | ~2.75% |

General Information

| | |
|---------------------|--|
| Storage | 25 to 40°C (Room Temperature) |
| Shelf Life | 60 Months |
| IMDG Identification | UN No.:1325 , IMCO Class No.:4.1 , Packing Group:III |
| HSN Code | |
| 1 Gms | 28030090 (GST 18%) |
| 250 Mg | 28030090 (GST 18%) |
| 50 Mg | 28030090 (GST 18%) |

Available Packages

50 Mg
250 Mg

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com