



Sisco Research
Laboratories
Pvt. Ltd.

Product Datasheet
and Specification

19093

Allyl Alcohol extrapure, 99%

Part A

CAS : 107-18-6

Molecular Formula : C3H6O

Molecular Weight : 58.08

Specifications

Appearance (Clarity)	Clear
Appearance (Colour)	Colourless
Appearance (Form)	Liquid
Assay (GC)	min. 99%
Density (g/ml) @ 20°C	0.850-0.852
Refractive Index (20°C)	1.411-1.413
Boiling Range	97-98°C

Other Information

Notes	(Forbidden by Air)
-------	--------------------

General Information

Storage	25 to 40°C (Room Temperature)
Shelf Life	60 Months
IMDG Identification	UN No.:1098 , IMCO Class No.:6.1(3) , Packing Group:I
HSN Code	
2500 ml	29052900 (GST 18%)
500 ml	29052900 (GST 18%)
Type of Packing	
2500 ml	Glass Bottle
500 ml	Glass Bottle

Available Packages

500 ml
2500 ml

Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srlchem.com