



19219

## Bradford Reagent for Proteins

Part B

### Specifications

Appearance (Colour)	Pale yellowish brown to brown to brownish green
Appearance (Form)	Solution
Suitable for protein determination	Passes test

### Other Information

Application	Bradford Reagent has been used to determine total protein concentration.
-------------	--

### General Information

Storage	2 to 8°C (Refrigerate)
Shelf Life	12 Months
IMDG Identification	UN No.:1805 , IMCO Class No.:8 , Packing Group:III
HSN Code	
100 ml	38229090 (GST 12%)
500 ml	38229090 (GST 12%)

### Available Packages

100 ml

---

500 ml

#### Disclaimer

The information represented here may/may not represent the entire product specification, application or protocol recommended by Sisco Research Laboratories Pvt. Ltd. (SRL). This information is for the user scientists or trading community as a guide in their applications. The company claims no liability for misuse resulting due to wrong usage of the information above. For actual batch related documents, mail us.

608-B, Satellite Gazebo, Andheri Ghatkopar Link Road, Chakala, Andheri (E), Mumbai - 400 099, Maharashtra, India. Telephone: +91-22-4268 5800, Email: info@srchem.com, website www.srchem.com